



Diabetes In-Home Care Kit

The Ultimate In-Home Care Starter Kit and In-Hospital Training Kit for Diabetic Pets

Diapetic™ Rescue provides everything a pet owner or veterinarian practice will need to address the post-diagnosis training and care of a diabetic pet. The ease-of-use can reduce the apprehension of caring for a diabetic pet and the tools contained in the kit provide a comfortable application for diabetes management. Kit testing components will last 80 days if testing daily. Refills can be ordered for replacement as needed. Insulin and Syringes sold separately via prescription through your Veterinarian. All components have been tested and endorsed by a veterinary teaching school. More information can be found at www.diapetic.com.

This kit is very useful for the newly diagnosed diabetic pet and their owner. You will learn to manage your pet's diabetes and be able to work with your veterinarian by learning to test and treat your beloved animal at home in a non-stressful environment. You and your pet will be happy to have help in a kit as you develop expertise in managing diabetes while saving time and money. The kit will be useful in managing your pet's illness for the long run in your home, while traveling or when training a pet-sitter.

PRODUCT FEATURES

- Convenient, Easy-to-Use Kit to Help You Navigate Diabetes
- Complete Training Instructions to Prepare You For Various Situation
- Consolidated Supply Management With Everything At Your Fingertips
- Tested with Proven Results Within a Renowned Veterinary School
- Endorsed by Endocrinologists and Veterinarians for Home or Training

www.diapetic.com
1-877-733-1870

COMPONENTS:

- Diabetic Pet Guidelines (Detailed laminated simple instructions for managing your pet's diabetes)
- Diabetic Pet ID Tag (Reflective red heart tag for both cats and dogs)
- Blood Glucose Monitor (Disposable monitor with 50 strips)
- Lancing Device (Gentle means to extract blood for testing)
- Keto-Diastix Urine Test (Recommended by all veterinarians with 50 strips)
- Honey Sticks (Use for low blood sugar)
- 12 Hour Ice Pack (Keeps insulin cool)
- Alcohol Swab (Sanitize injection area)
- Needle Disposal (Compact and safe needle disposal)
- Hand Sanitizer (Used by pet owner to ensure uncompromised testing results)
- Tissues (Used to wipe excess blood from ear & comfort back of ear during lancing)



A contribution will be made to UGA Arch Foundation towards improving treatment and diagnosis of diabetes in pets.

INSPIRE • CREATE • SIMPLIFY